

## Defibrillator checklist

Every year, ambulance services respond to about 60,000 cases of suspected cardiac arrests in the UK. Defibrillation within 3 – 5 min of collapse can produce survival rates as high as 50–70%. Every minute that defibrillation is delayed, a casualty's chances of survival falls by about 10%.\*

## We've put together a checklist to help you ensure your defibrillator is always ready for use.

Defibrillators self-test on a daily, weekly or monthly basis - so your device will signal, typically with a flashing light or audible alert, if there is a problem, for example, pads not correctly connected, or low battery.

Set up a regular inspection procedure for your defibrillator



**Record your inspections below** 



Record the expiration date of pads and **batteries** 



**Reorder batteries and pads** before they expire



Follow the guidance in your product manual when cleaning your defibrillator



| Defibrillator ready to use? | Date | Initial |
|-----------------------------|------|---------|
|                             |      |         |
|                             |      |         |
|                             |      |         |
|                             |      |         |



## Need new pads or batteries? Order now at sja.org.uk/shop or call 0344 770 4808†

2015 Resuscitation Council Guidelines. Resuscitation Council UK